

Lab Report on:

Determination of the Young's Modulus by the flexure of a beam

Experiment NU: 4

Course Title: Physics Lab I Course Code: PHY -102

Submitted by:

Md. Abu Raiyeen Reshad (ID: 20234103103) Intake: 52 (Section: 03), CSE

Submitted To:

Md. Sharif Ahmed, Lecturer, Department of EEE

Signature of Teacher

Submission Date:



Experiment NO: 04

Experiment Name! Determination of the Young's Modulus by the Plexure of a beam.

Apparatus:

- (1) Pin and Microscope
- 2. Meter Scale
- 3. Suitable weights
- 1. Slide caliper's

 - 5. Screw gauge 6. A long Metal beam
 - 7. Knife screw

Provided the distortion of a body is not too great it has been found that the amount of distortion is directly proportional to the magnitude of the forces producing the distortion. The fact is know as "Hooke's" law. If a wine of matural length i is structed on compressed a distance x by a force F, experiment reveals that


$$F = k\chi - - - (1)$$

where k is the constant whose value will depend on the material, the dimension of the wire and the units for measurement. The value of k is depend only on the material of the specimen and noton its dimension. It is called Young's Modulus of modulus of elascity for the modernial symbol y.

If a force F be applied normally to a crosssectional area A of the maderial in the front form of a wire, then F/A is called the tensile stress. x be the increase in length produced in an original length L as a result of this force, then x/1, is called the tensile strain.

$$\therefore Y = \frac{\text{densile stress}}{\text{densile strain}} = \frac{F/A}{\chi/L} = \frac{mg}{\pi n^{V}} \times \frac{L}{\chi} \frac{\text{dynes}}{\text{cm}^{2}}$$

If a rectangular beam of breadth band thickness d is supported near its two ends by two knife edges seperated by a distance 1, and if a mass m acting at a point of the beam equidistant from the knife edges producess a depression of then the Young's modulus is given by $Y = \frac{mg13}{4bd3x}$



$$Y = \frac{mgl^{3}}{4bd^{3}x}$$

$$= \frac{1.5 \times 8.8 \times (0.78)^{3}}{4 \times 0.038 \times (0.0038)^{3} \times 0.00403}$$

$$= 2.07 \times 10^{11} \text{Nm}^{2}$$
us by the flexure of a beam

Experiment No.: 04

Name of the Experiment: — — — — — — — — — — — Determination of the Young's Modulus by the flexure of a beam

Data Collection:

Length of the beam, $l = \frac{78.74}{1.000}$

Table 1: Data for load versus elongation

Additi	,	Readings for the elongation, x						Mean	Mean			
onal	Load increasing				Load decreasing				reading	depres		
Load	LSR	CSD	LC	CSR	Total	LSR	CSD	LC	CSR	Total	(cm)	sion
on	х	N	(cm)	$y = N \times$	Reading	х	N	(cm)	y = N	Reading		y_0
hange	(cm)			L.C.	=x+y	(cm)		1.7	$\times LC$	=x+y	,	(cm)
r (kg)				(cm)	cm)	1		(cm)	cm	. 1	
0	6	O		0	0	0.1	110		0.110	0.27	0.105	0.105
0.5	0.2	195		0.195	0.395	0.3	305	ς.	0.305	0.605	0.5	0.00
1.0	0.4		0.00%	0.438	0.838	0.6	550	200	0.550	1.15	0.999	0.495
1.5		613	0.6	0.613	1.313	0.8	670	The second	0.670	1.47	1.3915	
2.0	0.9	770		0.770	1.67	0.9	180	0	0.780	1.68	0 / /	
2,5	1.1	975		0.975	2.075	4.1	975		0.975	2.07.5	2.075	0.4
				<u>, </u>								

Table 2: Measure the breadth, (b) of beam

No. of obs.	Main scale reading (M.S.R) cm	Vernier scale divisions (V.S.D)	Vernier constant (V.C) cm	Vernier scale reading (V.S.R) = (V.S.D X V.C) cm	Total breadth b (cm)	Mean Breadth b (cm)
1	3.70	03	0.002	0.006	3.706	
2	3.80	21	0.002	0.092	38421	3.8226
3	3.90	lD	0.002	0.02	3.92	61

Table 3: Measure the depth, (d) of beam

X	No. of obs.	Main scale reading (M.S.R) cm	Vernier scale divisions (V.S.D)	Vernier constant (V.C) cm	Vernier scale reading (V.S.R) = (V.S.D X V.C) cm	Total depth d (cm)	Mean depth d (cm)
	l	0.3	40	0.002	0.08	0.380	
	2	0.3	43	0.002	0.086	0.386	0.380
1	3	0.3	37	0.002	0.074	0.379	\

From the sheet we get,

Length of the beam (1)= 78.74 cm = 0.78 m mean depression (x) = 0.4035 cm = 0.00403 m breadth of the beam (b)= 3.8226 cm = 0.038 m depth of the beam (d)= 0.380 = 0.0038 m given, weight of the beam (m)= 1.5 kg

we know,
$$Y = \frac{mgL^{3}}{4 bd^{3}\chi}$$

$$= \frac{1.5 \times 9.8 \times (0.78)^{3}}{4 \times (0.038) \times (0.0038)^{3} \times 0.00403}$$

Resulto

Young's modulus of the given material is 2.07 × 10" Nm-2

Precission:

- 1. First of all we have to find out the point equivalent reading of Main scale readin and circular scale, which we will count as o.
- 2. After increasing 0.5 kg in each steps, we have reset our circular, main scale. After reseting those scale we count every full rotation of our sircular scale.
 - 3. To find out the mean depression, we take the value which occurs friequently.
 - 1. For 2,5 kg the load increasing data and load decreasing data with be was some,

References:

- (i) Lab manual
- (ii) Links -> scrubd.com