<u>Name of the Experiment:</u>
Familiarization with Vernier Calipers and Screw Gauge by measuring various samples

 Table 1: Reading for Measurement of the length of the sample (Iron block) (using slide Calipers)

Dimension	No.	Main	Vernier	Vernier	Vernier scale	Total	Average	Volume of
to be	of	scale	scale	constant	reading	length/breadth/Th		the Iron
measured	obs	reading	divisions	(V.C)	(V.S.R) =	ickness, 1/b/h		block (V)
		(M.S.R)	(V.S.D)	(cm)	$(V.S.D \times V.C)$	(M.S.R + V.S.R)		cm <sup>3</sup>
		cm			cm)	(cm)		
Length of	1							
Iron block	2							
	3							
Breadth of Iron block	1							
HOH DIOCK	2							
	3							
Thickness of Iron	1							
block	2							
	3							

Table 2: Reading for Measurement of the radius of the sample (Lead Shot) (using Screw Gauge)

No. of obs.	Liner scale reading (L.S.R) cm	Circular scale divisions (C.S.D)	Least count (L.C) (cm)	Circular scale reading (C.S.R) = (C.S.D × L.C)	Total diameter (D) = (L.S.R + C.S.R) (cm)	Mean diameter D (cm)	Volume of the Lead Shot (V) (cm <sup>3</sup> )
1							
2							
3							

Name of the Experiment:
To determine the value of g, acceleration due to gravity, by means of a Compound Pendulum

Table 1: Data for time period

Starting direction	Hole no.	Distance from Top (cm)	Time for 10 oscillations (sec)	Mean Time, t (sec)	Mean Time Period, $T = \frac{t}{10}$ (sec)
	1	5	(i) (ii)		
	2	10	(i) (ii)		
	3	15	(i) (ii)		
	4	20	(i) (ii)		
Forward	5	25	(i) (ii)		
	6	30	(i) (ii)		
	7	35	(i) (ii)		
	8	40	(i) (ii)		
	9	45	(i) (ii)		
	1	55	(i) (ii)		
	2	60	(i) (ii)		
	3	65	(i) (ii)		
	4	70	(i) (ii)		
Reverse	5	75	(i) (ii)		
	6	80	(i) (ii)		
	7	85	(i) (ii)		
_	8	90	(i) (ii)		
	9	95	(i) (ii)		

Name of the Experiment:
To determine the spring constant and effective mass of a given spiral spring

 Table 1: Data of time period for different masses

No. of obs.	Loads m <sub>0</sub>	Extension 1	Times for 10 oscillation		Total Perio	Mean T (sec)	T <sup>2</sup>	
	(gm)	(cm)	$t_1$	$t_2$	$T_1 = \frac{t_1}{10}$	$T_2 = \frac{t_2}{10}$		
1								
2								
3								
4								
5								

# Name of the Experiment:

Determination of the Young's Modulus by the flexure of a beam

### **Data Collection:**

Length of the beam,  $l = \dots$  cm

 Table 1: Data for load versus elongation

Additi				Read	ings for the	elongatio	on, x				Mean	Mean
onal		I	Load incr	easing		Load decreasing					reading	depress
Load	L.S.R	C.S.D	L.C	C.S.R	Total	L.S.R	C.S.D	L.C	C.S.R	Total	(cm)	ion
on	(cm)		(cm)	=	Reading	(cm)		(cm)	=	Reading		$y_0$
hanger				$(C.S.D \times$	=(L.S.R				(C.S.D	=(L.S.R		(cm)
(kg)				L.C)	+				$\times$ L.C)	+		
				(cm)	C.S.R)				(cm)	C.S.R)		
					(cm)					(cm)		

**Table 2:** Measure the breadth, (b) of beam

No. of obs.	Main scale reading (M.S.R) (cm)	Vernier scale divisions (V.S.D)	Vernier constant (V.C) (cm)	Vernier scale reading (V.S.R) = (V.S.D × V.C) (cm)	Total breadth (b) = (MSR+VSR) (cm)	Mean Breadth b (cm)
1						
2						
3						

**Table 3:** Measure the depth, (d) of beam

No. of obs.	Main scale reading (M.S.R) (cm)	Vernier scale divisions (V.S.D)	Vernier constant (V.C) (cm)	Vernier scale reading (V.S.R) = (V.S.D × V.C) (cm)	$Total \\ depth \\ d = (MSR+VSR) \\ (cm)$	Mean depth d (cm)
1						
2						
3						

Name of the Experiment:
To determine the surface tension of water by capillary tube method

**Table 1:** Measurement of inside radius (r) of the tubes

	Readings in cm for the bore												
be no.	Deft side (X <sub>2</sub> )  VSR Total					Right side( X <sub>1</sub> )				Difference diameter	Radius r=D/2	Mean Radius r	
Tu	M.S.R (cm)	V.S.D	V.C (cm)	VSR = (V.S.D ×V.C) (cm)	Total Reading $(x_1) = (M.S.R+V.S.R)$ $(cm)$	M.S.R (cm)	V.S.D	V.C (cm)	VSR = (V.S.D ×V.C) (cm)	Total Reading $(x_2) =$ $(M.S.R+$ $V.S.R)$ $(cm)$	D=X <sub>2</sub> ~X <sub>1</sub> (cm) (cm)		
A													
В													

Table 2: Determination of the height of the column of water 'h'

No.		Reading	of the wa	ater meni	scus		Readi	ng at the	tip of po	inter	Height	Mean
of										h	Height	
obse											$=h_1\sim h_2$	h
rvati on	M.S.R (cm)	V.S.D	V.C (cm)	V.S.R = (V.S.D ×V.C) (cm)		M.S.R (cm)	V.S.D	V.C (cm)	V.S.R = (V.S.D ×V.C) (cm)	Total Reading (h2)= (M.S.R+ V.S.R) (cm)	(cm)	
1												
2												
3												
4												

### **Name of the Experiment:**

To determine the modulus of rigidity of a wire by Statical method using Barton's apparatus

**Table 1:** Diameter of the fly-wheel (heavy cylinder)

No. of obs.	Main scale reading (M.S.R) (cm)	Vernier scale divisions (V.S.D)	Vernier constant (V.C) (cm)	Vernier scale reading (V.S.R) = (V.S.D × V.C) (cm)	Total Diameter D = (MSR+VSR) (cm)	Mean Diameter d (cm)
1						
2						
3						

**Table 2:** Radius of the wire (using screw gauge)

No. of obs.	Liner scale reading (L.S.R) (cm)	Circular scale divisions (C.S.D)	Least count (L.C) (cm)	Circular scale reading (C.S.R) = (CSD × L.C)	Total diameter D = (LSR+CSR) (cm)	Mean diameter D (MSR+VSR) (cm)	Mean radius r=D/2 (cm)
1							
2							
3							

**Table 3:** Reading for load-twist graph

No. of	Load in each		Pointer read		Mean pointer reading		twist $\varphi 1^{\circ}$ )	
Obs.	hanger	Scale S <sub>1</sub>			$S_1$	$S_2$	n tw -φ1	
	(gm)	Load increasing	Load decreasing	Load increasing	Load decreasing	$(\varphi 1^{\circ})$	$(\varphi 2^{\circ})$	Mean (φ2°-φ
1								
2								
3								
4								
5								

# **Name of the Experiment:**

To determine the modulus of rigidity of a wire by method of oscillations (dynamic method)

### **Data collection:**

Length of the wire, L = .... cm

 Table 1: Readings for the Radius of the cylinder, R (using slide calipers)

No. of obs.	Main scale reading (M.S.R) (cm)	Vernier scale divisions (V.S.D)	Vernier constant (V.C) (cm)	Vernier scale reading (V.S.R) = (V.S.D × V.C) (cm)	Total diameter D = (MSR+VSR)	Mean diameter D (cm)	Mean radius R= D/2 (cm)
1							
2							
3							

**Table 2:** Radius of the wire, r (using screw gauge)

No. of obs.	Liner scale reading (L.S.R) (cm)	Circular scale divisions (C.S.D)	Least count (L.C) (cm)	Circular scale reading (C.S.R) = (CSD × L.C)	Total diameter D = (LSR+CSR) (cm)	Mean diameter D (cm)	Mean radius r=D/2 (cm)
1							
2							
3							

**Table 3:** Reading for the time period T.

No. of obs.	Time for 20 oscillations, t (see)	Period of oscillation $T = \frac{t}{20} (sec)$	Mean T. (sec)
1			
2			
3			

Experiment Name:

To determine the refractive index of the material of a prism

**Table 1:** Determination of angle of Prism

No. of obs.	Readings for image in the face AB of the prism				Readings for image in the face AC of the prism				Difference in readings at the two face $(\theta=M-N)$	Mean (θ) of the two venires	Angle of the prism $A=\theta/2$		
	M.S.R	V.S.D	V.C	Value of V.S.R	Total Reading N	M.S.R	V.S.D	V.C	Value of V.S.R	Total Reading M	Difference in the two face	Mean (θ ve	Angle of th
1													
2													
3													

Table 2: Determination of angle of minimum deviation

No. of obs.	Readings for the minimum deviation position					Readings for the direct position				Angle of minimum deviation $(\delta_m)$ = $(M-N)$	Mean $(\delta_m)$	
	M.S.R	V.S.D	V.C	Value of V.S.R	Total Reading N	M.S.R	V.S.D	V.C	Value of V.S.R	Total Reading M	Angle o deviation	Mea
1												
2												
3												

### **Name of the Experiment:**

To determine the value of g, acceleration due to gravity, by means of Kater's Pendulum

# **Data Collection:**

- **1.** Distance between  $K_1$  and  $CG(h_1)=.....cm$
- **2.** Distance between  $K_2$  and  $CG(h_2)=.....cm$
- **3.** Table for time period  $T_1$  (oscillation about  $K_1$ )

No. of	Number of	Time of Oscillation	Time Period	Mean Time
Observation	Oscillation, n	$t_1$ (sec)	$T_1 = \frac{t_1}{n}(\sec)$	$T_1(sec)$
1				
2				
3				
4				
5				

**4.** Table for time period  $T_2$  (oscillation about  $K_2$ )

No. of	Number of	Time of Oscillation	Time Period	Mean Time
Observation	Oscillation, n	$t_2$ (sec)	$T_2 = \frac{t_2}{n} \text{ (sec)}$	$T_2(sec)$
1				
2				
3				
4				
5				